

RotaCut[®]

Sphincterotomes - Single Use



Manufactured by:

Medi-Globe

Medi-Globe GmbH
Medi-Globe-Strasse 1-5
83101 Achenmühle, Germany
Phone: ++49 8032 973 379
Fax: ++49 8032 973 339
e-mail: sales@medi-globe.de
www.medi-globe.de

Distributed by:

Ver. 3.0 - 10/07

Medi-Globe
The Spirit of Care.™

RotaCut®

Sphincterotomes - Single Use

Features:

- 2 & 3 Lumen Versions
- Rotatable Tip
- Color Markings with Black Rings
- "Pre-curved" Versions
- Ergonomically Shaped Handle

Benefits:


- Contrast Medium Injection while Guide Wire is in Place
- Precise Positioning
- Macroscopically Controlled Cannulation Depth
- Easier Positioning to the Papilla
- User-Friendly Handling

Order Information:

3-Lumen Guide Wire Sphincterotomes, rotatable:

Catalog No.	Description	Outer Diameter	Cutting Wire Length	Guide Wire
GSP-31-20-020	Tapered Tip, Tip Length 5 mm	2.55/2.05 mm	20 mm	.035" (0.89 mm)
GSP-31-20-025	Tapered Tip, Tip Length 5 mm	2.55/2.05 mm	25 mm	.035" (0.89 mm)
GSP-31-20-030	Tapered Tip, Tip Length 5 mm	2.55/2.05 mm	30 mm	.035" (0.89 mm)

Tips:



2-Lumen Guide Wire Sphincterotomes, rotatable:

Catalog Nr.	Description	Outer Diameter	Cutting Wire Length	Guide Wire
GSP-21-17-020	Tapered Tip, Tip Length 5 mm	1.7 mm	20 mm	.021" (0.53 mm)
GSP-21-17-030	Tapered Tip, Tip Length 5 mm	1.7 mm	30 mm	.021" (0.53 mm)
GSP-22-17-020	Precut	1.7 mm	20 mm	.021" (0.53 mm)
GSP-22-17-030	Precut	1.7 mm	30 mm	.021" (0.53 mm)
GSP-21-20-020	Tapered Tip, Tip Length 5 mm	2.55/2.05 mm	20 mm	.035" (0.89 mm)
GSP-21-20-030	Tapered Tip, Tip Length 5 mm	2.55/2.05 mm	30 mm	.035" (0.89 mm)
GSP-24-20-020	Tapered Tip, Tip Length 10 mm	2.55/2.05 mm	20 mm	.035" (0.89 mm)
GSP-24-20-030	Tapered Tip, Tip Length 10 mm	2.55/2.05 mm	30 mm	.035" (0.89 mm)
GSP-23-20-030	Tapered Tip, Tip Length 20 mm	2.55/2.05 mm	30 mm	.035" (0.89 mm)
GSP-22-20-020	Precut	2.55/2.05 mm	20 mm	.035" (0.89 mm)
GSP-22-20-030	Precut	2.55/2.05 mm	30 mm	.035" (0.89 mm)
GSP-34-20-020	B II Sphincterotome	2.55/2.05 mm	Variable	.035" (0.89 mm)

Standard-Sphincterotomes, 1-Lumen, rotatable:

Catalog No.	Description	Outer Diameter	Cutting Wire Length	Tip Length
GSP-01-17-020	Standard Tip	1.7 mm	20 mm	5.0 mm
GSP-01-17-030	Standard Tip	1.7 mm	30 mm	5.0 mm
GSP-02-17-020	Precut	1.7 mm	20 mm	-
GSP-02-17-030	Precut	1.7 mm	30 mm	-
GSP-03-17-020	Precut with Short Tip	1.7 mm	20 mm	1.0 mm
GSP-03-17-030	Precut with Short Tip	1.7 mm	30 mm	1.0 mm
GSP-04-17-020	Tapered Tip	1.7 mm	20 mm	5.0 mm
GSP-04-17-030	Tapered Tip	1.7 mm	30 mm	5.0 mm
GSP-14-18-030	B II Sphincterotome	1.8 mm	Variable	-

Needle Knife Sphincterotomes, 1-Lumen, non-rotatable:

Catalog No.	Description	Outer Diameter
GSP-11-17-005	Needle Knife Sphincterotome, Needle Projection 0 - 5 mm	1.7 mm
GSP-11-17-010	Needle Knife Sphincterotome, Needle Projection 0 - 10 mm	1.7 mm
GSP-11-17-015	Needle Knife Sphincterotome, Needle Projection 0 - 15 mm	1.7 mm

Special-Sphincterotomes, non-rotatable:

Catalog No.	Description	Outer Diameter
GSP-12-18-004	Fistolotome; cutting spiral 4 mm; color markings	1.8 mm
GSP-13-17-015	Shark Fin Sphincterotome; cutting wire 15 mm; color markings	1.7 mm

Medi-Globe handles are directly compatible with the Olympus and electro-surgical generators. To ensure compatibility with other types of electro-surgical generators, please contact Medi-Globe GmbH or your distributor. We recommend the use of compatible Medi-Globe active cords. You will find the corresponding Medi-Globe active cords in the leaflet "Accessories".