

RotaCut[®] - M

Monofilament Sphincterotomes - Single Use

MONOFILAMENT CUTTING WIRE

- > Precise & accurate cut
- > Smooth-running sphincterotomy



3-Lumen Version

 **Medi-Globe[®]**
The Spirit of Care.™



RotaCut[®] - M

Monofilament Sphincterotomes - Single Use

Features:

Benefits

- Monofilament Cutting Wire -----> Precise & Smooth Incision of the Papillary Roof
- 3-Lumen Version -----> Contrast Medium Injection while Guide Wire is in Place
- Rotatable Cutting Wire -----> Precise Positioning
- Color Markings with Black Rings -----> Macroscopically Controlled Cannulation Depth
- "Pre-curved" Versions -----> Easier Positioning to the Papilla
- Ergonomically Shaped Handle -----> User-Friendly Handling

Order Information:

3-Lumen Guide Wire Sphincterotomes with monofilament cutting wire



Catalog No.	Description	Outer Diameter	Cutting Wire Length	Guide Wire
GSP-32-20-020	Tapered Tip, Tip Length 5 mm	2.55 / 2.05 mm	20 mm	.035" (0.89 mm)
GSP-32-20-025	Tapered Tip, Tip Length 5 mm	2.55 / 2.05 mm	25 mm	.035" (0.89 mm)
GSP-32-20-030	Tapered Tip, Tip Length 5 mm	2.55 / 2.05 mm	30 mm	.035" (0.89 mm)

Manufactured by:

 **Medi-Globe**

Medi-Globe GmbH
Medi-Globe-Strasse 1-5
83101 Achenmühle, Germany
Phone: ++49 8032 973 379
Fax: ++49 8032 973 339
e-mail: sales@medi-globe.de
www.medi-globe.de

Distributed by: