



ReBite®

➤ Detachable Cups

Optimal Cleaning & Enhanced Instrument Life

➤ Ground Distal Spiral (10 cm)

Easy Passage through Curved Endoscope Positions

Pulmonary

Oval Cups



ReBite® Biopsy Forceps with Fenestrated Cups

Catalog No.	Cup Shape	Outer Diameter	Working Length	Minimum Accessory Channel Diameter
401240	Oval	1.8 mm	120 cm	2.0 mm
406240	Oval	2.3 mm	120 cm	2.5 mm

Oval Cups with Spike



ReBite® Biopsy Forceps with Fenestrated Cups

Catalog No.	Cup Shape	Outer Diameter	Working Length	Minimum Accessory Channel Diameter
402240	Oval with Spike	1.8 mm	120 cm	2.0 mm
407240	Oval with Spike	2.3 mm	120 cm	2.5 mm

Serrated Oval Cups



ReBite® Biopsy Forceps with Fenestrated Cups

Catalog No.	Cup Shape	Outer Diameter	Working Length	Minimum Accessory Channel Diameter
405240	Serrated Oval Cups	1.8 mm	120 cm	2.0 mm
408240	Serrated Oval Cups	2.3 mm	120 cm	2.5 mm

Gastro & Colon

Oval Cups



ReBite® Biopsy Forceps with Fenestrated Cups

Catalog No.	Cup Shape	Outer Diameter	Working Length	Minimum Accessory Channel Diameter
401320	Oval	1.8 mm	160 cm	2.0 mm
401360	Oval	1.8 mm	180 cm	2.0 mm
401460	Oval	1.8 mm	230 cm	2.0 mm
406300	Oval	2.3 mm	150 cm	2.5 mm
406360	Oval	2.3 mm	180 cm	2.5 mm
406460	Oval	2.3 mm	230 cm	2.5 mm
406600	Oval	2.3 mm	300 cm	2.5 mm

Oval Cups with Spike



ReBite® Biopsy Forceps with Fenestrated Cups

Catalog No.	Cup Shape	Outer Diameter	Working Length	Minimum Accessory Channel Diameter
402320	Oval with Spike	1.8 mm	160 cm	2.0 mm
402360	Oval with Spike	1.8 mm	180 cm	2.0 mm
402460	Oval with Spike	1.8 mm	230 cm	2.0 mm
407320	Oval with Spike	2.3 mm	160 cm	2.5 mm
407360	Oval with Spike	2.3 mm	180 cm	2.5 mm
407460	Oval with Spike	2.3 mm	230 cm	2.5 mm

Serrated Oval Cups



ReBite® Biopsy Forceps with Fenestrated Cups

Catalog No.	Cup Shape	Outer Diameter	Working Length	Minimum Accessory Channel Diameter
405320	Serrated Oval Cups	1.8 mm	160 cm	2.0 mm
405360	Serrated Oval Cups	1.8 mm	180 cm	2.0 mm
405460	Serrated Oval Cups	1.8 mm	230 cm	2.0 mm
408320	Serrated Oval Cups	2.3 mm	160 cm	2.5 mm
408360	Serrated Oval Cups	2.3 mm	180 cm	2.5 mm
408460	Serrated Oval Cups	2.3 mm	230 cm	2.5 mm

Serrated Oval Cups with Spike



ReBite® Biopsy Forceps with Fenestrated Cups

Catalog No.	Cup Shape	Outer Diameter	Working Length	Minimum Accessory Channel Diameter
408360WN	Serrated Oval Cups with Spike	2.3 mm	180 cm	2.5 mm
408460WN	Serrated Oval Cups with Spike	2.3 mm	230 cm	2.5 mm

ReBite® Biopsy Forceps in different lengths as listed are available on request.

Special Versions

Oval Cups



ReBite® Mini Biopsy Forceps with Fenestrated Cups

Catalog No.	Cup Shape	Outer Diameter	Working Length	Minimum Accessory Channel Diameter
400220	Oval	1.0 mm	110 cm	1.2 mm
400240	Oval	1.0 mm	120 cm	1.2 mm
400360	Oval	1.0 mm	180 cm	1.2 mm
400460	Oval	1.0 mm	230 cm	1.2 mm

Oval Cups



Oval Cups with Spike



ReBite® Jumbo Biopsy Forceps with Extra Large Fenestrated Cups

Catalog No.	Cup Shape	Outer Diameter	Working Length	Minimum Accessory Channel Diameter
411360	Oval	3.4 mm	180 cm	3.6 mm
411460	Oval	3.4 mm	230 cm	3.6 mm
412360	Oval with Spike	3.4 mm	180 cm	3.6 mm
412460	Oval with Spike	3.4 mm	230 cm	3.6 mm

ReBite® - “Hot”

“Hot”-Biopsy Forceps – Reusable

ReBite®-“Hot”

Medi-Globe’s reusable “Hot” biopsy forceps is a two-piece instrument. One unit is comprised of the long catheter and forceps cups as referenced in our product catalog, (catalog no. 406XXXH & 407XXXH). The second required device component is the handle, (catalog no.100500) which is compatible with all Medi-Globe “Hot” biopsy forceps catheter units. This two-piece component system permits individual component replacement as necessary, without having to purchase a complete new instrument.

Medi-Globe’s “Hot” biopsy forceps handle is directly compatible with Olympus manufactured electro-surgical generators. To ensure compatibility with other types of electro-surgical generator units, please refer to the “Active Cords” section of this leaflet to select the appropriate active cord for your electro-surgical generator.

Oval Cups



Oval Cups with Spike



“Hot”-Biopsy Forceps (excluding handle), Oval Fenestrated Cups, Teflon Insulated Metal Spiral

Catalog No.	Cup Shape	Outer Diameter	Working Length	Minimum Accessory Channel Diameter
406360H	Oval	2.5 mm	180 cm	2.7 mm
406460H	Oval	2.5 mm	230 cm	2.7 mm
407360H	Oval Cups with Spike	2.5 mm	180 cm	2.7 mm
407460H	Oval Cups with Spike	2.5 mm	230 cm	2.7 mm

Handle with Connection to HF-Generators

Catalog No.	Description
100500	Handle for “Hot” - Biopsy Forceps