

Biopsy

Biopsie

i



 **Medi-Globe[®]**
The Spirit of Care.™

FlexiBite

Biopsy Forceps – Single Use

FlexiBite

- Sharp Standard & Serrated Cups
Precise & Clean Tissue Acquisition
- Spiked Versions
Controlled & Accurate Sampling from the Targeted Biopsy Site
- Versions with Coated Metal Spiral with Positioning Bands
Friction Minimized Passage through Working Channel
Convenient Evaluation of Forceps Position for a Controlled & Safe Catheter Advancement and Removal
- Hot-Biopsy Forceps
Safe Electro-coagulation and Excision of Small Polyps



 **Medi-Globe**[®]
The Spirit of Care.[™]

Cup Shapes



Biopsy Forceps Standard



FlexiBite – Biopsy Forceps with Metal Spiral, sterile packed

Catalog No.	Cup Shapes	Outer Diameter	Working Length	Minimum Working Channel Diameter
GBF-01-18-120	Oval Cups	1.8 mm	120 cm	2.0 mm
GBF-01-23-180	Oval Cups	2.3 mm	180 cm	2.5 mm
GBF-01-23-230	Oval Cups	2.3 mm	230 cm	2.5 mm
GBF-02-23-180	Oval Cups with Spike	2.3 mm	180 cm	2.5 mm
GBF-02-23-230	Oval Cups with Spike	2.3 mm	230 cm	2.5 mm



FlexiBite – Biopsy Forceps with Metal Spiral, sterile packed

Catalog No.	Cup Shapes	Outer Diameter	Working Length	Minimum Working Channel Diameter
GBF-03-23-230	Serrated Oval Cups	2.3 mm	230 cm	2.5 mm
GBF-04-23-230	Serrated Oval Cups with Spike	2.3 mm	230 cm	2.5 mm



Biopsy Forceps with Coated Metal Spiral

Oval Cups

Oval Cups
with Spike

FlexiBite – Biopsy Forceps with Coated Metal Spiral with Positioning Bands, sterile packed

Catalog No.	Cup Shapes	Outer Diameter	Working Length	Minimum Working Channel Diameter
GBF-21-23-180	Oval Cups	2.3 mm	180 cm	2.5 mm
GBF-21-23-230	Oval Cups	2.3 mm	230 cm	2.5 mm
GBF-22-23-180	Oval Cups with Spike	2.3 mm	180 cm	2.5 mm
GBF-22-23-230	Oval Cups with Spike	2.3 mm	230 cm	2.5 mm

Serrated Oval Cups

Serrated Oval Cups
with Spike

FlexiBite – Biopsy Forceps with Coated Metal Spiral with Positioning Bands, sterile packed

Catalog No.	Cup Shapes	Outer Diameter	Working Length	Minimum Working Channel Diameter
GBF-23-23-230	Serrated Oval Cups	2.3 mm	230 cm	2.5 mm
GBF-24-23-230	Serrated Oval Cups with Spike	2.3 mm	230 cm	2.5 mm



Positioning Bands on the Coating

“Hot”-Biopsy Forceps

Oval Cups



Monopolar “Hot”-Biopsy Forceps with Coated Metal Spiral, sterile packed



Catalog No.	Cup Shapes	Outer Diameter	Working Length	Minimum Working Channel Diameter
GBF-11-24-180	Oval Cups	2.4 mm	180 cm	2.8 mm
GBF-11-24-230	Oval Cups	2.4 mm	230 cm	2.8 mm

Medi-Globe “Hot”-Biopsy Forceps are directly compatible with Olympus electro-surgical generators. To verify compatibility with other electro-surgical generators please contact your Medi-Globe sales representative. Please also consult the Medi-Globe “Accessories” leaflet to see a complete list of active cords manufactured to fit most generators.

FlexiBite

Biopsy Forceps – Single Use

i



Distributed by:

Ver. 5.0 (2010-06)