




E.R.C.P. Catheters

- "Pre-Curved" Catheter with Memory Effect
Easier Positioning to the Papilla and Possible Manual Tip Manipulation for Individual Anatomy
- Color Markings
Macroscopically Controlled Cannulation Depth
- Metal Tip / Metal Rings
Maximum X-Ray Visualization

Tip Shapes



Compatible to .035" Guide Wire

| E.R.C.P. Catheter with Removable Nitinol Stylet | | |
|---|-------------------|-------------------------|
| Catalog No. | Tip Shape | Catheter Outer Diameter |
| 520420 | Standard Tip | 1.8 mm |
| 523420 | 10 mm Tapered Tip | 2.2 mm |
| 523421 | 5 mm Tapered Tip | 2.2 mm |
| 532420 | Bullet Metal Tip | 1.8 mm |
| 533420 | Ball Metal Tip | 1.8 mm |

Compatible to .021" Guide Wire

| E.R.C.P. Catheter with Removable Nitinol Stylet | | |
|---|-------------------|-------------------------|
| Catalog No. | Tip Shape | Catheter Outer Diameter |
| 522420 | 10 mm Tapered Tip | 1.8 mm |
| 522421 | 5 mm Tapered Tip | 1.8 mm |

Special Version

| Minor E.R.C.P. Catheter with Removable Nitinol Stylet | | |
|---|---------------------------|-------------------------|
| Catalog No. | Tip Shape | Catheter Outer Diameter |
| 534420 | Metal Cannula Length 1 mm | 1.8 mm |