




Drainage

Drainage


i



 **Medi-Globe[®]**
The Spirit of Care.™

Biliary & Pseudo-Cyst Drainage Stents

Single Use



Biliary & Pseudo-Cyst Drainage Stents

- Well Balanced Wall Diameter
Maximum Inner Lumen with Highest Form Stability & Kink Resistance
- Tapered & Round Ends
Atraumatic Implantation
- Pseudo-Cyst Double Pigtail Stents with Special Markers
Precise Placement via EUS Guidance
- Pre-Assembled Stent Introducer Sets
Convenient & Time-Saving Positioning of Biliary Stents



 **Medi-Globe®**
The Spirit of Care.™

Biliary & Pseudo-Cyst Drainage Stents

Medi-Globe offers the biliary stent as single component (catalog no. without "S") and additionally as integrated part of a set. Depending on the stent diameter and design sets consist of the following components:


7 French Stents and All Pigtail / Double Pigtail Stents

- Stent
- Pusher
- Teflon Coated Stainless Steel Guide Wire

Straight Stents in Diameter Larger than 7 French

- Stent
- Pusher
- Teflon Coated Stainless Steel Guide Wire
- Guiding Catheter


Medi-Globe manufactures stents in numerous sizes, which are not all listed in this leaflet. You are kindly requested to contact us to check the availability of the requested sizes.



Biliary Stent Set, Straight with Bottom Hole (Amsterdam Type)


Catalog No.	Stent		Guiding Catheter		Pusher		Guide Wire	
	French Size	Length between Flaps	French Size	Length	French Size	Length	Diameter	Length
GBS-01-07-005S	7.0	5 cm	N/A	N/A	7.0	220 cm	.035" (.089 mm)	450 cm
GBS-01-07-007S	7.0	7 cm	N/A	N/A	7.0	220 cm	.035" (.089 mm)	450 cm
GBS-01-07-009S	7.0	9 cm	N/A	N/A	7.0	220 cm	.035" (.089 mm)	450 cm
GBS-01-07-012S	7.0	12 cm	N/A	N/A	7.0	220 cm	.035" (.089 mm)	450 cm
GBS-01-07-015S	7.0	15 cm	N/A	N/A	7.0	220 cm	.035" (.089 mm)	450 cm
GBS-01-08-005S	8.5	5 cm	5.0	320 cm	8.5	160 cm	.035" (.089 mm)	450 cm
GBS-01-08-007S	8.5	7 cm	5.0	320 cm	8.5	160 cm	.035" (.089 mm)	450 cm
GBS-01-08-009S	8.5	9 cm	5.0	320 cm	8.5	160 cm	.035" (.089 mm)	450 cm
GBS-01-08-012S	8.5	12 cm	5.0	320 cm	8.5	160 cm	.035" (.089 mm)	450 cm
GBS-01-08-015S	8.5	15 cm	5.0	320 cm	8.5	160 cm	.035" (.089 mm)	450 cm
GBS-01-10-005S	10.0	5 cm	6.0	320 cm	10.0	160 cm	.035" (.089 mm)	450 cm
GBS-01-10-007S	10.0	7 cm	6.0	320 cm	10.0	160 cm	.035" (.089 mm)	450 cm
GBS-01-10-009S	10.0	9 cm	6.0	320 cm	10.0	160 cm	.035" (.089 mm)	450 cm
GBS-01-10-012S	10.0	12 cm	6.0	320 cm	10.0	160 cm	.035" (.089 mm)	450 cm
GBS-01-10-015S	10.0	15 cm	6.0	320 cm	10.0	160 cm	.035" (.089 mm)	450 cm
GBS-01-11-005S	11.5	5 cm	6.0	320 cm	11.5	160 cm	.035" (.089 mm)	450 cm
GBS-01-11-007S	11.5	7 cm	6.0	320 cm	11.5	160 cm	.035" (.089 mm)	450 cm
GBS-01-11-009S	11.5	9 cm	6.0	320 cm	11.5	160 cm	.035" (.089 mm)	450 cm
GBS-01-11-012S	11.5	12 cm	6.0	320 cm	11.5	160 cm	.035" (.089 mm)	450 cm
GBS-01-11-015S	11.5	15 cm	6.0	320 cm	11.5	160 cm	.035" (.089 mm)	450 cm





Biliary Stent Set, Straight with Side Ports

Catalog No.	Stent		Guiding Catheter		Pusher		Guide Wire	
	French Size	Length between Flaps	French Size	Length	French Size	Length	Diameter	Length
GBS-02-07-004S	7.0	4 cm	N/A	N/A	7.0	220 cm	.035" (.089 mm)	450 cm
GBS-02-07-006S	7.0	6 cm	N/A	N/A	7.0	220 cm	.035" (.089 mm)	450 cm
GBS-02-07-008S	7.0	8 cm	N/A	N/A	7.0	220 cm	.035" (.089 mm)	450 cm
GBS-02-07-011S	7.0	11 cm	N/A	N/A	7.0	220 cm	.035" (.089 mm)	450 cm
GBS-02-07-016S	7.0	16 cm	N/A	N/A	7.0	220 cm	.035" (.089 mm)	450 cm
GBS-02-08-004S	8.5	4 cm	5.0	320 cm	8.5	160 cm	.035" (.089 mm)	450 cm
GBS-02-08-006S	8.5	6 cm	5.0	320 cm	8.5	160 cm	.035" (.089 mm)	450 cm
GBS-02-08-008S	8.5	8 cm	5.0	320 cm	8.5	160 cm	.035" (.089 mm)	450 cm
GBS-02-08-011S	8.5	11 cm	5.0	320 cm	8.5	160 cm	.035" (.089 mm)	450 cm
GBS-02-08-016S	8.5	16 cm	5.0	320 cm	8.5	160 cm	.035" (.089 mm)	450 cm
GBS-02-10-004S	10.0	4 cm	6.0	320 cm	10.0	160 cm	.035" (.089 mm)	450 cm
GBS-02-10-006S	10.0	6 cm	6.0	320 cm	10.0	160 cm	.035" (.089 mm)	450 cm
GBS-02-10-008S	10.0	8 cm	6.0	320 cm	10.0	160 cm	.035" (.089 mm)	450 cm
GBS-02-10-011S	10.0	11 cm	6.0	320 cm	10.0	160 cm	.035" (.089 mm)	450 cm
GBS-02-10-016S	10.0	16 cm	6.0	320 cm	10.0	160 cm	.035" (.089 mm)	450 cm
GBS-02-11-004S	11.5	4 cm	6.0	320 cm	11.5	160 cm	.035" (.089 mm)	450 cm
GBS-02-11-006S	11.5	6 cm	6.0	320 cm	11.5	160 cm	.035" (.089 mm)	450 cm
GBS-02-11-008S	11.5	8 cm	6.0	320 cm	11.5	160 cm	.035" (.089 mm)	450 cm
GBS-02-11-011S	11.5	11 cm	6.0	320 cm	11.5	160 cm	.035" (.089 mm)	450 cm
GBS-02-11-016S	11.5	16 cm	6.0	320 cm	11.5	160 cm	.035" (.089 mm)	450 cm




Biliary Stent Set, Maxi-Flow

Catalog No.	Stent		Guiding Catheter		Pusher		Guide Wire	
	French Size	Length between Flaps	French Size	Length	French Size	Length	Diameter	Length
GBS-03-07-005S	7.0	5 cm	N/A	N/A	7.0	220 cm	.035" (.089 mm)	450 cm
GBS-03-07-007S	7.0	7 cm	N/A	N/A	7.0	220 cm	.035" (.089 mm)	450 cm
GBS-03-07-009S	7.0	9 cm	N/A	N/A	7.0	220 cm	.035" (.089 mm)	450 cm
GBS-03-07-012S	7.0	12 cm	N/A	N/A	7.0	220 cm	.035" (.089 mm)	450 cm
GBS-03-07-015S	7.0	15 cm	N/A	N/A	7.0	220 cm	.035" (.089 mm)	450 cm
GBS-03-08-005S	8.5	5 cm	5.0	320 cm	8.5	160 cm	.035" (.089 mm)	450 cm
GBS-03-08-007S	8.5	7 cm	5.0	320 cm	8.5	160 cm	.035" (.089 mm)	450 cm
GBS-03-08-009S	8.5	9 cm	5.0	320 cm	8.5	160 cm	.035" (.089 mm)	450 cm
GBS-03-08-012S	8.5	12 cm	5.0	320 cm	8.5	160 cm	.035" (.089 mm)	450 cm
GBS-03-08-015S	8.5	15 cm	5.0	320 cm	8.5	160 cm	.035" (.089 mm)	450 cm


Biliary & Pseudo-Cyst Drainage Stents

Single Use




Biliary Stent Set, Maxi-Flow


Catalog No.	Stent		Guiding Catheter		Pusher		Guide Wire	
	French Size	Length between Flaps	French Size	Length	French Size	Length	Diameter	Length
GBS-03-10-005S	10.0	5 cm	6.0	320 cm	10.0	160 cm	.035" (.089 mm)	450 cm
GBS-03-10-007S	10.0	7 cm	6.0	320 cm	10.0	160 cm	.035" (.089 mm)	450 cm
GBS-03-10-009S	10.0	9 cm	6.0	320 cm	10.0	160 cm	.035" (.089 mm)	450 cm
GBS-03-10-012S	10.0	12 cm	6.0	320 cm	10.0	160 cm	.035" (.089 mm)	450 cm
GBS-03-10-015S	10.0	15 cm	6.0	320 cm	10.0	160 cm	.035" (.089 mm)	450 cm
GBS-03-11-005S	11.5	5 cm	6.0	320 cm	11.5	160 cm	.035" (.089 mm)	450 cm
GBS-03-11-007S	11.5	7 cm	6.0	320 cm	11.5	160 cm	.035" (.089 mm)	450 cm
GBS-03-11-009S	11.5	9 cm	6.0	320 cm	11.5	160 cm	.035" (.089 mm)	450 cm
GBS-03-11-012S	11.5	12 cm	6.0	320 cm	11.5	160 cm	.035" (.089 mm)	450 cm
GBS-03-11-015S	11.5	15 cm	6.0	320 cm	11.5	160 cm	.035" (.089 mm)	450 cm



Biliary Stent Set, Multi-Flap


Catalog No.	Stent		Guiding Catheter		Pusher		Guide Wire	
	French Size	Length between Flaps	French Size	Length	French Size	Length	Diameter	Length
GBS-04-07-005S	7.0	5 cm	N/A	N/A	7.0	220 cm	.035" (.089 mm)	450 cm
GBS-04-07-007S	7.0	7 cm	N/A	N/A	7.0	220 cm	.035" (.089 mm)	450 cm
GBS-04-07-009S	7.0	9 cm	N/A	N/A	7.0	220 cm	.035" (.089 mm)	450 cm
GBS-04-07-012S	7.0	12 cm	N/A	N/A	7.0	220 cm	.035" (.089 mm)	450 cm
GBS-04-07-015S	7.0	15 cm	N/A	N/A	7.0	220 cm	.035" (.089 mm)	450 cm
GBS-04-08-005S	8.5	5 cm	5.0	320 cm	8.5	160 cm	.035" (.089 mm)	450 cm
GBS-04-08-007S	8.5	7 cm	5.0	320 cm	8.5	160 cm	.035" (.089 mm)	450 cm
GBS-04-08-009S	8.5	9 cm	5.0	320 cm	8.5	160 cm	.035" (.089 mm)	450 cm
GBS-04-08-012S	8.5	12 cm	5.0	320 cm	8.5	160 cm	.035" (.089 mm)	450 cm
GBS-04-08-015S	8.5	15 cm	5.0	320 cm	8.5	160 cm	.035" (.089 mm)	450 cm
GBS-04-10-005S	10.0	5 cm	6.0	320 cm	10.0	160 cm	.035" (.089 mm)	450 cm
GBS-04-10-007S	10.0	7 cm	6.0	320 cm	10.0	160 cm	.035" (.089 mm)	450 cm
GBS-04-10-009S	10.0	9 cm	6.0	320 cm	10.0	160 cm	.035" (.089 mm)	450 cm
GBS-04-10-012S	10.0	12 cm	6.0	320 cm	10.0	160 cm	.035" (.089 mm)	450 cm
GBS-04-10-015S	10.0	15 cm	6.0	320 cm	10.0	160 cm	.035" (.089 mm)	450 cm
GBS-04-11-005S	11.5	5 cm	6.0	320 cm	11.5	160 cm	.035" (.089 mm)	450 cm
GBS-04-11-007S	11.5	7 cm	6.0	320 cm	11.5	160 cm	.035" (.089 mm)	450 cm
GBS-04-11-009S	11.5	9 cm	6.0	320 cm	11.5	160 cm	.035" (.089 mm)	450 cm
GBS-04-11-012S	11.5	12 cm	6.0	320 cm	11.5	160 cm	.035" (.089 mm)	450 cm
GBS-04-11-015S	11.5	15 cm	6.0	320 cm	11.5	160 cm	.035" (.089 mm)	450 cm





Biliary Stent Set, Curved

Catalog No.	Stent		Guiding Catheter		Pusher		Guide Wire	
	French Size	Length between Flaps	French Size	Length	French Size	Length	Diameter	Length
GBS-05-07-005S	7.0	5 cm	N/A	N/A	7.0	220 cm	.035" (.089 mm)	450 cm
GBS-05-07-007S	7.0	7 cm	N/A	N/A	7.0	220 cm	.035" (.089 mm)	450 cm
GBS-05-07-009S	7.0	9 cm	N/A	N/A	7.0	220 cm	.035" (.089 mm)	450 cm
GBS-05-07-012S	7.0	12 cm	N/A	N/A	7.0	220 cm	.035" (.089 mm)	450 cm
GBS-05-07-015S	7.0	15 cm	N/A	N/A	7.0	220 cm	.035" (.089 mm)	450 cm
GBS-05-08-005S	8.5	5 cm	5.0	320 cm	8.5	160 cm	.035" (.089 mm)	450 cm
GBS-05-08-007S	8.5	7 cm	5.0	320 cm	8.5	160 cm	.035" (.089 mm)	450 cm
GBS-05-08-009S	8.5	9 cm	5.0	320 cm	8.5	160 cm	.035" (.089 mm)	450 cm
GBS-05-08-012S	8.5	12 cm	5.0	320 cm	8.5	160 cm	.035" (.089 mm)	450 cm
GBS-05-08-015S	8.5	15 cm	5.0	320 cm	8.5	160 cm	.035" (.089 mm)	450 cm
GBS-05-10-005S	10.0	5 cm	6.0	320 cm	10.0	160 cm	.035" (.089 mm)	450 cm
GBS-05-10-007S	10.0	7 cm	6.0	320 cm	10.0	160 cm	.035" (.089 mm)	450 cm
GBS-05-10-009S	10.0	9 cm	6.0	320 cm	10.0	160 cm	.035" (.089 mm)	450 cm
GBS-05-10-012S	10.0	12 cm	6.0	320 cm	10.0	160 cm	.035" (.089 mm)	450 cm
GBS-05-10-015S	10.0	15 cm	6.0	320 cm	10.0	160 cm	.035" (.089 mm)	450 cm
GBS-05-11-005S	11.5	5 cm	6.0	320 cm	11.5	160 cm	.035" (.089 mm)	450 cm
GBS-05-11-007S	11.5	7 cm	6.0	320 cm	11.5	160 cm	.035" (.089 mm)	450 cm
GBS-05-11-009S	11.5	9 cm	6.0	320 cm	11.5	160 cm	.035" (.089 mm)	450 cm
GBS-05-11-012S	11.5	12 cm	6.0	320 cm	11.5	160 cm	.035" (.089 mm)	450 cm
GBS-05-11-015S	11.5	15 cm	6.0	320 cm	11.5	160 cm	.035" (.089 mm)	450 cm




Biliary Stent Set, Single Pigtail

Catalog No.	Stent		Pusher		Guiding Wire	
	French Size	Length between Pigtail	French Size	Length	Diameter	Length
GBS-06-05-015S	5.0	15 cm	5.0	220 cm	.035" (0.89 mm)	450 cm
GBS-06-07-015S	7.0	15 cm	7.0	220 cm	.035" (0.89 mm)	450 cm
GBS-06-08-015S	8.5	15 cm	8.5	220 cm	.035" (0.89 mm)	450 cm
GBS-06-10-015S	10.0	15 cm	10.0	220 cm	.052" (1.32 mm)	400 cm


Biliary & Pseudo-Cyst Drainage Stents

Single Use




Biliary Stent Set, Double Pigtail

Catalog No.	Stent		Pusher		Guiding Wire	
	French Size	Length between Pigtails	French Size	Length	Diameter	Length
GBS-07-05-006S	5.0	6 cm	5.0	220 cm	.035" (0.89 mm)	450 cm
GBS-07-05-009S	5.0	9 cm	5.0	220 cm	.035" (0.89 mm)	450 cm
GBS-07-05-012S	5.0	12 cm	5.0	220 cm	.035" (0.89 mm)	450 cm
GBS-07-07-006S	7.0	6 cm	7.0	220 cm	.035" (0.89 mm)	450 cm
GBS-07-07-009S	7.0	9 cm	7.0	220 cm	.035" (0.89 mm)	450 cm
GBS-07-07-012S	7.0	12 cm	7.0	220 cm	.035" (0.89 mm)	450 cm
GBS-07-08-006S	8.5	6 cm	8.5	220 cm	.035" (0.89 mm)	450 cm
GBS-07-08-009S	8.5	9 cm	8.5	220 cm	.035" (0.89 mm)	450 cm
GBS-07-08-012S	8.5	12 cm	8.5	220 cm	.035" (0.89 mm)	450 cm
GBS-07-10-006S	10.0	6 cm	10.0	220 cm	.052" (1.32 mm)	400 cm
GBS-07-10-009S	10.0	9 cm	10.0	220 cm	.052" (1.32 mm)	400 cm
GBS-07-10-012S	10.0	12 cm	10.0	220 cm	.052" (1.32 mm)	400 cm



Pseudo-Cyst Double Pigtail Stents with Position Markers

Catalog No.	Stent	
	French Size	Length between Pigtails
GBS-11-08-002	8.5	2 cm
GBS-11-08-005	8.5	5 cm
GBS-11-08-007	8.5	7 cm
GBS-11-08-009	8.5	9 cm
GBS-11-10-002	10.0	2 cm
GBS-11-10-005	10.0	5 cm
GBS-11-10-007	10.0	7 cm
GBS-11-10-009	10.0	9 cm



Pre-assembled Stent-Introducer Sets

Pre-assembled Stent Introducer Set consists of a guiding catheter and pusher in one device. The set is a convenient option to single components for the placement of straight and curved biliary stents (catalog model: GBS-01-XX-XXX – GBS-05-XX-XXX), which require a guiding catheter. The set has an extra port to inject contrast medium during the entire procedure and while a guide wire is in place. Radiopaque metal markers on the guiding catheter and the pusher allow a precise positioning under X-ray view.



Pre-assembled Stent-Introducer Set with 1 Metal Ring on the Guiding Catheter

Catalog No.	Pusher			Guiding Catheter	
	French Size	Length between Flaps	Color Concept	French Size	Length
GBS-31-08-180	8.5	180 cm	orange	5.0	220 cm
GBS-31-10-180	10.0	180 cm	red	6.0	220 cm
GBS-31-11-180	11.5	180 cm	blue	6.0	220 cm



Pre-assembled Stent-Introducer Set with 4 Metal Rings in a Distance of 5 cm on the Guiding Catheter

Catalog No.	Pusher			Guiding Catheter	
	French Size	Length between Flaps	Color Concept	French Size	Length
GBS-32-08-180	8.5	180 cm	orange	5.0	220 cm
GBS-32-10-180	10.0	180 cm	red	6.0	220 cm
GBS-32-11-180	11.5	180 cm	blue	6.0	220 cm



Stent-Introducer Components

Single Use

Guiding Catheters

Guiding Catheters are available with one and with four metal rings. The metal markers on the 4 ring version are in a distance of 5 cm and allow a pre-measurement of the stricture length and ease the selection of the correct stent.



Guiding Catheter with 1 Metal Ring

Catalog No.	French Size	Length
GBS-41-05-220	5.0	220 cm
GBS-41-05-320	5.0	320 cm
GBS-41-06-220	6.0	220 cm
GBS-41-06-320	6.0	320 cm



Guiding Catheter with 4 Metal Rings

Catalog No.	French Size	Length
GBS-43-06-220	6.0	220 cm
GBS-43-06-320	6.0	320 cm

Pushers

Pushers follow for each diameter a defined color concept to differentiate the related french sizes.



Pusher with 1 Metal Ring

Catalog No.	French Size	Color Concept	Length
GBS-42-05-220	5.0	grey	220 cm
GBS-42-07-160	7.0	yellow	160 cm
GBS-42-07-220	7.0	yellow	220 cm
GBS-42-08-160	8.5	orange	160 cm
GBS-42-08-180	8.5	orange	180 cm
GBS-42-08-220	8.5	orange	220 cm
GBS-42-10-160	10.0	red	160 cm
GBS-42-10-180	10.0	red	180 cm
GBS-42-10-220	10.0	red	220 cm
GBS-42-11-160	11.5	blue	160 cm
GBS-42-11-180	11.5	blue	180 cm
GBS-42-11-220	11.5	blue	220 cm

Distributed by:

Manufactured by:

Medi-Globe GmbH | Medi-Globe-Straße 1-5 | 83101 Achenmühle, Germany
 ☎ +49 80 32 / 973- 379 | 📠 +49 80 32 / 973- 339 | ✉ sales@medi-globe.de

www.medi-globe.de

Ver. 5.0 (2010-10)

