

Pancreatic Stents

Single Use


Pancreatic Stents

- Well Balanced Wall Diameter
*Maximum Inner Lumen with Highest Form
Stability & Kink Resistance*
- Tapered & Round Ends
Atraumatic Implantation
- Radiopaque Material
Exact Control of Placement



 **Medi-Globe**[®]
The Spirit of Care.[™]

i



Pancreatic Stents


Medi-Globe offers the pancreatic stents as single component (catalog no. without "S") and additionally as integrated part of a set. Depending on the stent diameter and shape sets include the following components:

7 French Stents and All Pigtail / Double Pigtail Stents


- Stent
- Pusher
- Teflon Coated Stainless Steel Guide Wire

Straight Stents in Diameter Larger than 7 French

- Stent
- Pusher
- Teflon Coated Stainless Steel Guide Wire
- Guiding Catheter




Catalog No.	Stent		Pusher		Guide Wire	
	French Size	Length between Pigtail	French Size	Length	Diameter	Length
GBS-08-07-006S	7.0	6 cm	7.0	220 cm	.035" (0.89 mm)	450 cm
GBS-08-07-009S	7.0	9 cm	7.0	220 cm	.035" (0.89 mm)	450 cm
GBS-08-07-012S	7.0	12 cm	7.0	220 cm	.035" (0.89 mm)	450 cm
GBS-08-08-006S	8.5	6 cm	8.5	220 cm	.035" (0.89 mm)	450 cm
GBS-08-08-009S	8.5	9 cm	8.5	220 cm	.035" (0.89 mm)	450 cm
GBS-08-08-012S	8.5	12 cm	8.5	220 cm	.035" (0.89 mm)	450 cm
GBS-08-10-006S	10.0	6 cm	10.0	220 cm	.052" (1.32 mm)	400 cm
GBS-08-10-009S	10.0	9 cm	10.0	220 cm	.052" (1.32 mm)	400 cm
GBS-08-10-012S	10.0	12 cm	10.0	220 cm	.052" (1.32 mm)	400 cm
GBS-08-11-006S	11.5	6 cm	11.5	220 cm	.052" (1.32 mm)	400 cm
GBS-08-11-009S	11.5	9 cm	11.5	220 cm	.052" (1.32 mm)	400 cm
GBS-08-11-012S	11.5	12 cm	11.5	220 cm	.052" (1.32 mm)	400 cm



Pancreatic Stent Straight with Radial Holes		
Catalog No.	French Size	Length between Flap
GBS-10-07-005	7.0	5 cm
GBS-10-07-007	7.0	7 cm
GBS-10-07-008	7.0	8 cm
GBS-10-07-009	7.0	9 cm
GBS-10-07-011	7.0	11 cm
GBS-10-07-016	7.0	16 cm

Pancreatic Stent Straight with Radial Holes		
Catalog No.	French size	Length between Flap
GBS-10-08-005	8.5	5 cm
GBS-10-08-007	8.5	7 cm
GBS-10-08-008	8.5	8 cm
GBS-10-08-009	8.5	9 cm
GBS-10-08-011	8.5	11 cm
GBS-10-08-016	8.5	16 cm





Pancreatic Stent Set S - Shape

Catalog No.	Stent		Guiding Catheter		Pusher		Guide Wire	
	French Size	Length	French Size	Length	French Size	Length	Diameter	Length
GBS-09-05-005S	5.0	5 cm	N/A	N/A	5.0	220 cm	.035" (.089 mm)	450 cm
GBS-09-05-007S	5.0	7 cm	N/A	N/A	5.0	220 cm	.035" (.089 mm)	450 cm
GBS-09-05-009S	5.0	9 cm	N/A	N/A	5.0	220 cm	.035" (.089 mm)	450 cm
GBS-09-05-012S	5.0	12 cm	N/A	N/A	5.0	220 cm	.035" (.089 mm)	450 cm
GBS-09-05-015S	5.0	15 cm	N/A	N/A	5.0	220 cm	.035" (.089 mm)	450 cm
GBS-09-07-005S	7.0	5 cm	N/A	N/A	7.0	220 cm	.035" (.089 mm)	450 cm
GBS-09-07-007S	7.0	7 cm	N/A	N/A	7.0	220 cm	.035" (.089 mm)	450 cm
GBS-09-07-009S	7.0	9 cm	N/A	N/A	7.0	220 cm	.035" (.089 mm)	450 cm
GBS-09-07-012S	7.0	12 cm	N/A	N/A	7.0	220 cm	.035" (.089 mm)	450 cm
GBS-09-07-015S	7.0	15 cm	N/A	N/A	7.0	220 cm	.035" (.089 mm)	450 cm
GBS-09-10-005S	10.0	5 cm	6.0	320 cm	10.0	220 cm	.035" (.089 mm)	450 cm
GBS-09-10-007S	10.0	7 cm	6.0	320 cm	10.0	220 cm	.035" (.089 mm)	450 cm
GBS-09-10-009S	10.0	9 cm	6.0	320 cm	10.0	220 cm	.035" (.089 mm)	450 cm
GBS-09-10-012S	10.0	12 cm	6.0	320 cm	10.0	220 cm	.035" (.089 mm)	450 cm
GBS-09-10-015S	10.0	15 cm	6.0	320 cm	10.0	220 cm	.035" (.089 mm)	450 cm

Guiding Catheters

Guiding Catheters are available with one and with four metal rings. The metal markers on the 4 ring version are in a distance of 5 cm and allow a pre-measurement of the stricture length and ease the selection of the correct stent.



Guiding Catheter with 1 Metal Ring

Catalog No.	French Size	Length
GBS-41-05-220	5.0	220 cm
GBS-41-05-320	5.0	320 cm
GBS-41-06-220	6.0	220 cm
GBS-41-06-320	6.0	320 cm



Guiding Catheter with 4 Metal Rings

Catalog No.	French Size	Length
GBS-43-06-220	6.0	220 cm
GBS-43-06-320	6.0	320 cm


Pushers

Pushers follow for each diameter a defined color concept to differentiate the related French sizes.



Pusher with 1 Metal Ring



Catalog No.	French Size	Color Concept	Length
GBS-42-05-220	5.0	grey	220 cm
GBS-42-07-160	7.0	yellow	160 cm
GBS-42-07-220	7.0	yellow	220 cm
GBS-42-08-160	8.5	orange	160 cm
GBS-42-08-180	8.5	orange	180 cm
GBS-42-08-220	8.5	orange	220 cm
GBS-42-10-160	10.0	red	160 cm
GBS-42-10-180	10.0	red	180 cm
GBS-42-10-220	10.0	red	220 cm
GBS-42-11-160	11.5	blue	160 cm
GBS-42-11-180	11.5	blue	180 cm
GBS-42-11-220	11.5	blue	220 cm



Pancreatic Stents

Single Use

i



Distributed by:

Ver. 3.0 (2010-10)

www.medi-globe.de

Manufactured by:

Medi-Globe GmbH | Medi-Globe-Straße 1-5 | 83101 Achenmühle, Germany
☎ +49 80 32/973-379 | 📠 +49 80 32/973-339 | ✉ sales@medi-globe.de