

Aspiration Needles

for Bronchoscopy & Gastroscopy - Single Use

Aspiration Needles

- Coaxial Luer-Lock for Vacuum Syringe
Easy Aspiration of Cells
- Safety Clip and Lock Mechanism
Secure Needle Advancement and Projection
- 15° Beveled Needle with Side Hole
Atraumatic Piercing Capability
Maximum Yield of Aspirates



Aspiration Needles

for Bronchoscopy & Gastroscopy - Single Use

Bronchoscopy

Aspiration Needle – with Metal Spiral and including 10 ml Vacuum Syringe with Stopcock

| Catalog No. | Needle Diameter | Needle Length | Metal Spiral Outer Diameter | Working Length | Minimum Accessory Channel |
|---------------|-------------------|---------------|-----------------------------|----------------|---------------------------|
| GAN-01-19-100 | 21 Gauge / 0.8 mm | 12 mm | 1.9 mm | 100 cm | 2.1 mm |

Gastroscopy

Aspiration Needle – with Plastic Tube and including 10 ml Vacuum Syringe with Stopcock

| Catalog No. | Needle Diameter | Needle Length | Sheath Outer Diameter | Working Length | Minimum Accessory Channel |
|---------------|-------------------|---------------|-----------------------|----------------|---------------------------|
| GAN-02-19-180 | 21 Gauge / 0.8 mm | 12 mm | 1.9 mm | 180 cm | 2.1 mm |

Distributed by:

Manufactured by: