

Fibrin Needles

Single Use

Fibrin Needles

- Double Lumen Coaxiale „Needle in Needle“-System
Precise Confluence of Fibrinogen/Thrombin at Intended Position
- Removable Safety Clip
Secure Needle within the Catheter
- 15° Beveled Needle
Effective Needle Puncture
- Easy Grip Handle
Easy Handling



Fibrin Needles

Single Use

Fibrin Needles

Fibrin Needle				
Catalog No.	Needle Diameter	Needle Length	Catheter Outer Diameter	Working Length
GFN-01-04-180	1.2 mm	4.0 mm	2.5 mm	180 cm
GFN-01-04-230	1.2 mm	4.0 mm	2.5 mm	230 cm
GFN-01-06-180	1.2 mm	6.0 mm	2.5 mm	180 cm
GFN-01-06-230	1.2 mm	6.0 mm	2.5 mm	230 cm
GFN-01-08-180	1.2 mm	8.0 mm	2.5 mm	180 cm
GFN-01-08-230	1.2 mm	8.0 mm	2.5 mm	230 cm

Distributed by:

Manufactured by: