

# Polyp & Foreign Body Retrievers

Reusable



## Polyp & Foreign Body Retrievers

- Wide Opening of Graspers  
*Large Variety of Applications*
- Detachable Distal Tips (Forceps)  
*Optimal Cleaning & Enhanced Instrument Life*
- Ground Distal Spiral (10 cm)  
*Easy Passage through Curved Endoscope Positions*

4-Prong Graspers with Eyelets



### Polyp Retrievers

Catalog No.	Grasper Shape	Outer Diameter	Working Length
805360	4-Prong Grasper with Eyelets	2.5 mm	180 cm
805460	4-Prong Grasper with Eyelets	2.5 mm	230 cm

4-Prong Graspers with Hooks



### Foreign Body Retrievers

Catalog No.	Grasper Shape	Outer Diameter	Working Length
804360	4-Prong Grasper with Hooks	2.5 mm	180 cm
804460	4-Prong Grasper with Hooks	2.5 mm	230 cm

Forked Jaws 2:2



Forked Jaws 1:2



### Foreign Body Retrievers

Catalog No.	Grasper Shape	Outer Diameter	Working Length
807360	Forked Jaws 2:2	2.5 mm	180 cm
807460	Forked Jaws 2:2	2.5 mm	230 cm
812360	Forked Jaws 1:2	2.5 mm	180 cm
812460	Forked Jaws 1:2	2.5 mm	230 cm

Rat Tooth



Alligator



### Foreign Body Forceps with Detachable Tip

Catalog No.	Tip Shape	Outer Diameter	Working Length
816360	Rat Tooth	1.8 mm	180 cm
816460	Rat Tooth	1.8 mm	230 cm
817360	Rat Tooth	2.5 mm	180 cm
817460	Rat Tooth	2.5 mm	230 cm

### Foreign Body Forceps with Detachable Tip

Catalog No.	Tip Shape	Outer Diameter	Working Length
820360	Alligator	1.8 mm	180 cm
820460	Alligator	1.8 mm	230 cm
821360	Alligator	2.5 mm	180 cm
821460	Alligator	2.5 mm	230 cm