

Foreign Body Retrieval

Fremdkörperentnahme

i



 **Medi-Globe**[®]
The Spirit of Care.[™]



Polyp & Foreign Body Retrievers


Single Use

Polyp & Foreign Body Retrievers

- Strong & Wide Opening Graspers
Large Variety of Applications
- Ergonomically Shaped Handle
User-Friendly and Safe



 **Medi-Globe**[®]
The Spirit of Care.[™]



Polyp & Foreign Body Retrievers

Single Use

i

Polyp & Foreign Body Retrievers

Foreign Body Retriever



Grasper with Hooks

Polyp Retriever



Grasper with Eyelets

Foreign Body Retriever



Grasper with Forked Jaws (2:2)



Grasper with Forked Jaws (1:2)

Foreign Body Basket



Foreign Body Retriever – Grasper with Hooks

Catalog No.	Grasper Shape	Catheter Outer Diameter	Working Length
GFB-01-25-180	4-Prong Grasper	2.5 mm	180 cm
GFB-01-25-230	4-Prong Grasper	2.5 mm	230 cm

Polyp Retriever – Grasper with Eyelets

Catalog No.	Grasper Shape	Catheter Outer Diameter	Working Length
GFB-06-25-180	3-Prong Grasper	2.5 mm	180 cm
GFB-06-25-230	3-Prong Grasper	2.5 mm	230 cm
GFB-02-25-180	4-Prong Grasper	2.5 mm	180 cm
GFB-02-25-230	4-Prong Grasper	2.5 mm	230 cm

Foreign Body Retriever – Grasper with Forked Jaws

Catalog No.	Grasper Shape	Catheter Outer Diameter	Working Length
GFB-03-25-180	Forked Jaws 2:2	2.5 mm	180 cm
GFB-03-25-230	Forked Jaws 2:2	2.5 mm	230 cm
GFB-04-25-180	Forked Jaws 1:2	2.5 mm	180 cm
GFB-04-25-230	Forked Jaws 1:2	2.5 mm	230 cm

Foreign Body Basket

Catalog No.	Grasper Shape	Basket Diameter	Catheter Outer Diameter	Working Length
GFB-07-25-425	4-Wire Basket	25 mm	2.5 mm	230 cm
GFB-07-25-625	6-Wire Basket	25 mm	2.5 mm	230 cm

Ver. 3.0 (2010-08)

