



Stone Extraction Baskets

Medi-Globe's stone extraction basket sets offer a certain compatibility with the StoneBuster – mechanical lithotripsy system. Sets with the catalog model nos. 60XXXX are fully detachable and in emergency case the Teflon tube can be replaced against an emergency spiral (catalog no. GML-90-34-080) and connected to StoneBuster handle (catalog no. GML-99-01-000) to exercise an emergency lithotripsy for a basket removal out of the patient.

Combi - Sets with the catalog model nos. 615XXX are for use as stone extraction basket and depending on the situation also as lithotripsy basket. The Teflon tube can be replaced by a standard lithotripsy metal spiral (catalog no. GML-90-26-180) and in connection with StoneBuster handle (catalog no. GML-99-01-000) a through the endoscope mechanical lithotripsy can be performed.

Stone Extraction Basket Set – Consists of Basket, Outer Sheath (Working Length 180 cm), Stepless Adjustable Sealing Cap and Instrument Holder

Catalog no.	Description	Basket Diameter	Outer Sheath Diameter
606420	4-Wire Basket	20 mm	1.9 mm
606425	4-Wire Basket	25 mm	1.9 mm
606430	4-Wire Basket	30 mm	1.9 mm
608620	6-Wire Basket	20 mm	1.9 mm
608625	6-Wire Basket	25 mm	1.9 mm
608630	6-Wire Basket	30 mm	1.9 mm
605420	4-Wire Basket	20 mm	2.5 mm
605425	4-Wire Basket	25 mm	2.5 mm
605430	4-Wire Basket	30 mm	2.5 mm
607620	6-Wire Basket	20 mm	2.5 mm
607625	6-Wire Basket	25 mm	2.5 mm
607630	6-Wire Basket	30 mm	2.5 mm

Combi-Set Stone Extraction / Lithotripsy Basket – Consists of Basket, Outer Sheath (Working Length 280 cm), Stepless Adjustable Sealing Cap and Instrument Holder

Catalog no.	Description	Basket Diameter	Outer Sheath Diameter
615420	4-Wire Basket	20 mm	2.5 mm
615425	4-Wire Basket	25 mm	2.5 mm
615430	4-Wire Basket	30 mm	2.5 mm
615620	6-Wire Basket	20 mm	2.5 mm
615625	6-Wire Basket	25 mm	2.5 mm

All sets are equipped with an instrument holder and are ready for use. For convenience reasons Medi-Globe offers optionally a sturdy handle (catalog no. 100600) with thumb ring and finger slide.

Handle for Stone Extraction Baskets (Optionally)

Catalog No.	Description
100600	Handle for Stone Extraction Baskets