



RotaSnare®

Medi-Globe's reusable polypectomy snare is a three-piece instrument. One complete unit is comprised of

- a) Snare Wire
- b) Insertion Tube
- c) Handle for Polypectomy Snares (catalog no. 100700)

This component system permits individual component replacement as necessary, without having to purchase a complete new instrument.


Medi-Globe's handle for polypectomy snares is directly compatible with Olympus manufactured electro-surgical generators. To ensure compatibility with other types of electro-surgical generator units, please refer to the "Active Cords" section of this leaflet to select the appropriate active cord for your electro-surgical generator.

↗ 360° Rotating Versions
Easier Polyp Snaring

↗ Multifilament Loop Wire
Safe Cut and Long Lasting Loop

↗ Complete Component Detachability
Economical Component Replacement

Rotating Polypectomy Snares



RotaSnare® – with Multifilament Loop (Snare Wire without Handle and Insertion Tube)

Catalog No.	Loop Shape	Loop Diameter	Working Length	Compatible Insertion Tube Catalog No.
713360	Oval	15 mm	180 cm	713363
713460	Oval	15 mm	230 cm	713463
715360	Oval	25 mm	180 cm	713363
715460	Oval	25 mm	230 cm	713463
716360	Oval, with Insulated Tip	25 mm	180 cm	713363
716460	Oval, with Insulated Tip	25 mm	230 cm	713463
717360	Oval	35 mm	180 cm	713363
717460	Oval	35 mm	230 cm	713463
718360	Oval, with Insulated Tip	35 mm	180 cm	713363
718460	Oval, with Insulated Tip	35 mm	230 cm	713463



RotaSnare® – with Multifilament Loop (Snare Wire without Handle and Insertion Tube)

Catalog No.	Loop Shape	Loop Diameter	Working Length	Compatible Insertion Tube Catalog No.
721360	Hexagonal	25 mm	180 cm	713363
721460	Hexagonal	25 mm	230 cm	713463
722360	Hexagonal, with Insulated Tip	25 mm	180 cm	713363
722460	Hexagonal, with Insulated Tip	25 mm	230 cm	713463

Standard Polypectomy Snares, non-rotatable

Onion



Standard with Multifilament Loop (Snare Wire without Handle and Insertion Tube)

Catalog No.	Loop Shape	Loop Diameter	Working Length	Compatible Insertion Tube Catalog No.
723360	Onion	50 mm	180 cm	713363
723460	Onion	50 mm	230 cm	713463

Oval Monofilament



Standard with Monofilament Loop (Snare Wire without Handle and Insertion Tube)

Catalog No.	Loop Shape	Loop Diameter	Working Length	Compatible Insertion Tube Catalog No.
726360	Oval	25 mm	180 cm	713363
726460	Oval	25 mm	230 cm	713463

Hexagonal



Standard with Multifilament Loop (Snare Wire without Handle and Insertion Tube)

Catalog No.	Loop Shape	Loop Diameter	Working Length	Compatible Insertion Tube Catalog No.
720360	Hexagonal	25 mm	180 cm	713363
720460	Hexagonal	25 mm	230 cm	713463
725360	Hexagonal	35 mm	180 cm	713363
725460	Hexagonal	35 mm	230 cm	713463

Asymmetric



Standard with Multifilament Loop (Snare Wire without Handle and Insertion Tube)

Catalog No.	Loop Shape	Loop Diameter	Working Length	Compatible Insertion Tube Catalog No.
703360	Asymmetric	15 mm	180 cm	703363
703460	Asymmetric	15 mm	230 cm	703463
706360	Asymmetric	25 mm	180 cm	703363
706460	Asymmetric	25 mm	230 cm	703463

Insertion Tubes with Sheath Diameter 2.5 mm

Catalog No.	Compatible to Snare Wire catalog nos.
713363	713360, 715360, 716360, 717360, 718360, 720360, 721360, 722360, 723360, 725360, 726360
713463	713460, 715460, 716460, 717460, 718460, 720460, 721460, 722460, 723460, 725460, 726460
703363	703360, 706360
703463	703460, 706430

Handle with Connection to HF - Generators

Catalog No.	Description
100700	Handle for Polypectomy Snares