

Cleaning Brushes

Single Use

Cleaning Brushes

- Double-Sided Brushes
Efficient Cleaning of Residues
- Channel-Friendly Plastic Coated Brush Tip
Avoids Accessory Channel Scratch
- Colored Catheter
Quick Identification of Different Brush Sizes
- Stackable Dispenser Box
Hygienic and User-Friendly Storage



Cleaning Brushes

Single Use

Cleaning Brushes

Cleaning Brush – Double Sided, Individually Packed in a Dispenser Box

Catalog No.	Brush Diameter	Catheter Length	Accessory Channel Compatibility	Unit	Product
GAC-02-03-230	3 mm	230 cm	2.0 mm	Box of 100	
GAC-02-05-230	5 mm	230 cm	2.8 mm	Box of 100	
GAC-03-05-011	Brush 1: 5 mm / Brush 2: 11 mm	230 cm	2.8 - 6.0 mm	Box of 100	

Valve Cleaning Brush – Double Sided

Catalog No.	Brush Diameter	Unit	Product
GAC-03-05-012	Brush 1: 5 mm / Brush 2: 12 mm	Box of 25	

Dispenser Box



Distributed by: