

EasyPass Guide Wires

Single Use

EasyPass Guide Wires

- Super Flexible, Hydrophilic EasyPass Tip
Providing Excellent Accessibility & Guide Wire Passage through especially Challenging Anatomy
- High Performance Radiopaque Tip
Precise Guide Wire Control & Maneuverability under Fluoroscopic Imaging
- Two-Color Teflon Coated Nitinol Wire
Kink Resistant
May Be Left in Place During Sphincterotomy



 **Medi-Globe**[®]
The Spirit of Care.[™]

EasyPass Guide Wires

Single Use

EasyPass Guide Wires

EasyPass Guide Wires are featured with a super flexible 8 cm long highly radiopaque and hydrophilic tip. This EasyPass guide wire has outstanding performance characteristics and provides enhanced maneuverability through difficult to access anatomical structures.

The Teflon coated Nitinol guide wire may be left in place during sphincterotomy. A shorter length guide wire is offered for use in combination with the available short wire systems.

EasyPass Guide Wire – Two-Color Teflon Coated Nitinol Wire with 8 cm Super Flexible, Hydrophilic & Radiopaque Tip

Catalog No.	Description	Diameter	Length
GGW-10-35-260	“Super Flexible” Tip	.035”	260 cm
GGW-10-35-450	“Super Flexible” Tip	.035”	450 cm