



Luna

Ambulatory Urodynamics

Luna





The Luna is a very easy to use Ambulatory Urodynamic recorder.

Using Bluetooth® technique the Luna offers an on-line view of recordings.

The Luna is affordable and weighs less than 200 grams. Ambulatory Urodynamics has never been so easy.



The Luna can be hooked up with a wireless flowmeter using Bluetooth®



MMS can offer T-DOC Air-charged and solid state catheters from French sizes 3 up to 9.

- **Ambulatory Investigations**

Wireless recording of 3 Pressure channels (Bladder, Urethra and Abdominal), Uroflow, EMG and Automated Leak Point Pressures through Conductivity.

- **On-line View of Recordings**

Wireless Bluetooth® technology makes it possible to view data while recording. High quality signal recording can be guaranteed before patient leaves investigation room.

- **Complete Range of Catheters**

MMS has a complete range of single use T-DOC and reusable microtip catheters for use with Luna. Even small sizes like 3 or 4 French for children are available.

- **Lightweight Housing, Easy Carrying**

The Luna weighs less than 200 grams and can easily be carried around by little children making ambulatory urodynamic recording more reliable than ever.

- **Fast Presentation of Investigation Results**

More than 24 hour of investigation data can be downloaded, analyzed and presented in less than 2 minutes. All Detrusor-flow plots are standard available. Ambulatory Urodynamics has never been so easy.

MMS
Medical Measurement Systems

P.O. Box 580

7500 AN Enschede

The Netherlands

T + 31 - 53 - 480 37 00

F + 31 - 53 - 480 37 01

E info@mmsinternational.com

I www.mmsinternational.com