




Flowmaster

Wireless Uroflowmeter

Flowmaster



The Flowmaster is a wireless digital flowmeter designed for practical, everyday flow study investigations, where ease of use, simplicity and cost effectiveness are top priorities.



Flowmaster Software running on standard Microsoft Windows® based PC

- **Automatic Investigation and Analysis**
Investigation automatically starts when patient begins voiding.
- **Wireless data transmission**
No more messy cables, easy to move around.
- **Patient Database**
Secure storage of investigation data in Flowmaster Database. Possibility to store investigation data from max. 7 Flowmasters (used simultaneously) in one central Database.
- **Automatic Backup**
Possibility to store investigation data on the hospital server for automatic backup.
- **Digital Flowsensor with Artefact Detection**
Weighted transducer type flowsensor with intelligent software which reduces artefacts.
- **Compact**
Flowmaster is a small lightweight and portable system. Available with height adjustable flowstand or micturation chair.


Medical Measurement Systems

P.O. Box 580
7500 AN Enschede
The Netherlands
T + 31 - 53 - 480 37 00
F + 31 - 53 - 480 37 01
E info@mmsinternational.com
I www.mmsinternational.com